

Isle of Wight NHS

Continence Appliance Formulary

Community

Isle of Wight NHS Continence Service

Updated Sept 2019

CONTENTS

<i>Page Number</i>	<i>Content</i>
2.	Introduction
3.	Urethral Catheters
4.	Suprapubic Catheters
5.	Catheter Valves
6.	Leg & Belly Bags
7.	Night Bags
8.	Catheter Fixation Catheter Bag Sleeves
9.	Maintenance Solution
10.	Catheter Trays
11.	Urinary Sheaths Leg Bag for Sheaths
12.	Night Bags for Sheaths Skin Preparation for Sheath Use Sheath Fixation Devices Alternatives to Sheaths
13.	Intermittent Catheters
14.	Intermittent Catheters continued Meatal Dilators Anal Irrigation
15.	Ordering

Introduction

The purpose of this formulary is to present a range of products which should be used for catheter/sheath users. It has been written with both clinical quality and cost effectiveness in mind.

It has been RAG'd with the following colours and rationale:

Colour	Rationale
	First line use.
	Second line, should be used with caution and in some cases only be used once all green options have been tried and there is a justified clinical reason
	Only use with continence service advice
	Continence Service Use only

There is also a guide in each section for the recommended amounts to order; initial & repeat ordering (this excludes emergency call outs).

The price is included as a guide and is from the June 2019 drug tariff; therefore it is subject to change.

Catheters

Urinary & Suprapubic Catheters

Clinical patient assessment is essential before deciding to catheterise, due to the high risk of infection associated with catheter use. Always consider alternative ways of management before resorting to an indwelling catheter.

Except in exceptional circumstances, no more than a 14Ch should be used for female urethral catheterisation and no more than a 16Ch standard catheter for a male. ***Female length catheters should never be used urethrally in the male.***

*Lubricants are not listed; instillagel is only to be used, as per trust policy.

Urethral Catheters

Manufacturer	Product	Length / Size / Balloon size	Price (£)	Order Code
L.IN.C Medical Systems Ltd	All Silicone Catheter with UniBal type balloon with sterile water filled syringe for balloon inflation & empty syringe for balloon deflation	Standard 5ml / 12ch	5.57	08501205G
		5ml / 14ch	5.57	08501405G
	All Silicone Open Tipped Catheter with UniBal type balloon with sterile water filled syringe for balloon inflation & empty syringe for balloon deflation	Female 5ml / 12ch	5.57	085012051G
		5ml / 14ch	5.57	085014051G
Unomedical Ltd	Careline 100% Silicone Foley Catheter 2 way with pre-filled syringe	Standard 5ml / 12ch	10.95	08451205W-08451405W
		5ml / 14ch	10.95	08451205W-08451405W
	Female 5ml / 12ch	10.95	08471205W-08471405W	
		5ml / 14ch	10.95	08471205W-08471405W
Teleflex	Rusch Brilliant Aquaflate All silicone with sterile water filled syringe for balloon inflation & empty syringe for balloon deflation	Male 10ml / 12ch	6.11	DA310112 - DA310124
		10ml / 14ch	6.11	DA310112 - DA310124
	Female 10ml / 12ch	6.11	DA210112 - DA210124	
		10ml / 14ch	6.11	DA210112 - DA210124
L.IN.C Medical Systems Ltd	All Silicone Catheter with glycerine filled syringe for balloon inflation & empty syringe for balloon deflation	Tiemann Tip 5ml / 12ch 5ml / 14ch	8.57 8.57	1011-12-1011-18 1011-12-1011-18
Bard Ltd	Bard Biocath Hydrogel Coated	Male 10ml / 12ch	8.50	D2265
		10ml / 14ch	8.50	D2265
		30ml / 12ch	8.61	D2266
		30ml / 14ch	8.61	D2266
		30ml / 16ch	8.61	D2266
		Female 10ml / 12ch	8.61	D2269
		10ml / 14ch	8.61	D2269
		Unomedical Ltd	Careline 100% Silicone Foley Catheter 2 way with pre-filled syringe	Tiemann Tip 10 ml / 12ch
10ml / 14ch	6.50			UM41161210PS-UM41162410PS
10ml / 16ch	6.50			UM41161210PS-UM41162410PS
Bard Ltd	Bard Bardex I.C. with silver alloy coating and a pre-filled syringe of sterile watersterile	Standard 10ml / 12ch	9.92	D236512S-D236522S
		10ml / 14ch	9.92	D236512S-D236522S
		Female 10ml / 12ch	9.92	D236912S-D236916S
		10ml / 14ch	9.92	D236912S-D236916S

Recommended Order Amounts**Initial Order:** 3**Subsequent Orders:** when down to 1 catheter, reorder 3

Suprapubic Catheters

Manufacturer	Product	Length / Size / Balloon size	Price (£)	Order Code	
L.I.N.C Medical Systems Ltd	All Silicone Catheter with UniBal type balloon with sterile water filled syringe for balloon inflation & empty syringe for balloon deflation	Standard			
		5ml / 12ch	5.57	08501205G	
		5ml / 14ch	5.57	08501405G	
		10ml / 16ch	5.57	08501610G	
		Female			
		5ml / 12ch	5.57	085012051G	
	5ml / 14ch	5.57	085014051G		
	10ml / 16ch	5.57	085016101G		
	All Silicone Open Tipped Catheter with UniBal type balloon with sterile water filled syringe for balloon inflation & empty syringe for balloon deflation	Standard			
		5ml / 12ch	10.95	08451205W-08451405W	
5ml / 14ch		10.95	08451205W-08451405W		
10ml / 16ch		10.95	08451610W-08452410W		
Female					
5ml / 12ch		10.95	08471205W-08471405W		
5ml / 14ch	10.95	08471205W-08471405W			
10ml / 16ch	10.95	08471610W-08472010W			
Teleflex	Rusch Brilliant Aquafate All silicone with sterile water filled syringe for balloon inflation & empty syringe for balloon deflation	Male			
		10ml / 12ch	6.11	DA310112 - DA310124	
		10ml / 14ch	6.11	DA310112 - DA310124	
		10ml / 16ch	6.11	DA310112 - DA310124	
		Female			
		10ml / 12ch	6.11	DA210112 - DA210124	
10ml / 14ch	6.11	DA210112 - DA210124			
Unomedical Ltd	Careline 100% Silicone Foley Catheter 2 way with pre-filled syringe	Male			
		10ml / 12ch	5.00	UM41151210PS-UM41152410PS	
		10ml / 14ch	5.00	UM41151210PS-UM41152410PS	
		10ml / 16ch	5.00	UM41151210PS-UM41152410PS	
		Female			
		10ml / 12ch	5.00	UM41171210PS-UM41172410PS	
10ml / 14ch	5.00	UM41171210PS-UM41172410PS			
Bard Ltd	Bard Bardex I.C. with silver alloy coating and a pre-filled syringe	Standard			
		10ml / 12ch	9.92	D236512S-D236522S	
		10ml / 14ch	9.92	D236512S-D236522S	
		10ml / 16ch	9.92	D236512S-D236522S	
		Female			
		10ml / 12ch	9.92	D236912S-D236916S	
10ml / 14ch	9.92	D236912S-D236916S			
Bard Ltd	Bard Biocath Hydrogel Coated	Male			
		10ml / 12ch	8.50	D2265	
		10ml / 14ch	8.50	D2265	
		10ml / 16ch	8.50	D2265	
		30ml / 12ch	8.61	D2266	
		30ml / 14ch	8.61	D2266	
		30ml / 16ch	8.61	D2266	
		Female			
		10ml / 12ch	8.61	D2269	
		10ml / 14ch	8.61	D2269	
		10ml / 16ch	8.61	D2269	

All available in size 18-22ch, if patient has this size refer to continence service

Recommended Order Amounts

Initial Order: 3

Subsequent Orders: when down to 1 catheter, reorder 3

Drainage Systems

Maintenance of a closed system, using a sterile bag or catheter valve, is important in reducing infection rates.

Bags and valves connected directly to the catheter should be changed every 5 to 7 days, except the Rusch belly bag, which is changed every 28 days.

Catheter valves may be preferred by the patient and help to maintain bladder wall elasticity.

Catheter Valves

Are a comfortable and discreet alternative to leg bags and allow a “normal” bladder filling/emptying cycle.

They are not suitable where bladder distension is contra indicated (e.g. post-radical prostatectomy, renal impairment and ureteric reflux)

Valves are also not suitable for patients with small bladder capacity, absence of bladder sensation, cognitive impairment, or insufficient manual dexterity to operate valve.

Catheter Valves

Manufacturer	Product	Pack Size	Price (£)	Order Code
L.I.N.C Medical Systems Ltd	Care-Flo Catheter Valve	5	7.59	CF1
Bard Ltd	Flip-Flo Catheter Valve	5	13.48	BFF5

Recommended Oder Amounts

Initial Order: 2 packs of 5 (10)

Subsequent Orders: when down to 2 valves, reorder 2 packs

Leg & Belly Bags

Bags positioned more than 30cm below the bladder may cause the mucosa to be sucked into the catheter eyes leading to bladder trauma and problems with drainage.

An option for patients for whom leg bags and catheter valves are not suitable is the ‘Belly Bag’, although it is not recommended where the catheter blocks more often than every 4 weeks, or for patients who need a larger capacity than 1 litre overnight (as it does not connect to a night drainage bag).

Leg & Belly Bags

Manufacturer	Product	Tube / Tap	Pack Size	Price (£)	Order Code	
L.IN.C Medical Systems Ltd	LINC-Flo sterile leg drainage bag with overnight connector and 10 pairs of elastic velcro straps. 350ml	Direct / T-tap	10	20.72	LM350SD-T	
		10cm / T-tap	10	20.72	LM350MD-T	
	Linc-flo sterile leg bags with overnight connector and 10 pairs of elastic velcro straps. 500mls	Direct / T-tap	10	22.31	LM350SD-T	
		10cm / T-tap	10	22.31	LM500MD-T	
		30cm / T-tap	10	22.31	LM500LD-T	
		Adjustable / T-tap	10	22.31	LM500AD-T	
		Direct / Lever	10	22.27	LM500SD-L	
		10cm / Lever	10	22.27	LM500MD-L	
		30cm / Lever	10	22.27	LM500LD-L	
		Adjustable / Lever	10	22.01	LM500AD-T	
	Linc-flo sterile leg bags with overnight connector and 10 pairs of elastic velcro straps. 750mls	Direct / T-tap	10	22.28	LM750SD-T	
		10cm / T-tap	10	22.28	LM750MD-T	
		30cm / T-tap	10	22.28	LM750LD-T	
		30cm / T-tap	10	22.28	LM750AD-T	
		Adjustable / Lever	10	22.28	LM750AD-L	
	Teleflex	Rusch Leg Bag 350ml	10cm / slide	10	21.90	850361-000350
10cm / lever			10	21.90	850363-000350	
30cm / slide			10	21.90	850362-000350	
30cm / lever			10	21.90	850364-000350	
Rusch Leg Bag 500ml		10cm / slide	10	23.94	850361-000500	
		10cm / lever	10	23.94	850363-000500	
		30cm / slide	10	23.94	850362-000500	
		30cm / lever	10	23.94	850364-000500	
Rusch Leg Bag 750ml		10cm / slide	10	23.94	850363-000750	
		10cm / lever	10	23.94	850364-000750	
Rusch belly bag		Includes belt & can be used for 28 days	1	10.78	B1000	
Rusch belly bag with extended drainage tubing		Includes belt & can be used for 28 days	1	10.78	B1000CT	
Rusch belly bag with sample port		Includes belt & can be used for 28 days	1	10.78	B1000P	
Coloplast		Simpla Plus leg bag with two-step safety lock tap anti-kink tubing	350ml / 6cm	10	30.78	21561
			350ml / 25cm	10	30.78	21562
			500ml / 6cm	10	30.78	21571
			500ml / 25cm	10	30.78	21572
			500ml / 50cm	10	30.11	21583
	750ml / 6cm		10	30.78	21591	
	750ml / 25cm		10	30.78	21592	
	Simpla Profile leg bag with lever lock tap and anti-kink tubing		350ml / 6cm	10	30.34	21563
		350ml / 25cm	10	30.34	21564	
		500ml / 6cm	10	30.86	21573	
		500ml / 25cm	10	30.86	21574	
		500ml / 50cm	10	30.86	21584	
		750ml / 6cm	10	31.31	21593	
		750ml / 25cm	10	31.31	21594	
	Simpla Trident T1 with slide action tap	350ml / short	10	30.57	370802	
		500ml / long	10	30.80	370817	
		500ml / short	10	30.80	370807	
		750ml / long	10	30.95	370819	
		750ml / short	10	30.95	370809	
		750ml / adjustable long	10	30.95	370904	
Simpla Trident T2 with lever action tap	350ml / short	10	30.93	376137		
	500ml / short	10	31.27	376138		
	500ml / long	10	31.27	376139		
	750ml / short	10	31.42	376140		
	750ml / long	10	31.42	376142		
Unomedical Ltd	Careline Urine Drainage bag with lever tap	500ml / short	10	25.48	46-05-SVCG	
		500ml / long	10	25.48	46-06-LVCG	

Recommended Order Amounts

Initial Order: 1 box of 10

Subsequent Orders: reorder 1 box when 2 bags are left

Night Bags.

Patients who wear a leg bag or catheter valve, and who need additional capacity at night, should have a ‘**single use**’ night bag fitted in tandem.

Drainable (or re-useable) night bags are suitable for catheterised patients who are nursed in bed and do not need a leg bag. They are also appropriate for urinary sheath wearers.

The night bag should be supported on a stand, or hanger, not more than 30cm below bladder level. These are complimentary items and can be requested from any company when placing an order.

Night Bags

Manufacturer	Product	Capacity / Tap	Pack Size	Price (£)	Order Code
L.IN.C Medical Systems Ltd	LINC-Flo non-sterile overnight drainage bag with 100cm inlet tube	2L / twist off tube outlet	10	2.30	LM2LNS
	LINC-Flo sterile overnight drainage bag with 100cm inlet tube	2L / twist off tube outlet	10	2.65	LM2LS
Teleflex	Non-sterile night bag with single use	2 Litre / drainable tap	30	9.55	850422-000000
	Rusch night drainage bag	2 Litre / slide 2 Litre / lever	10 10	10.78 10.78	850420-000000 850421-000000
Unomedical Ltd	Careline Easi MT drainage bag with 90cm inlet tube	2L / twist off tube outlet	10	3.40	47-60-LBH

Recommended Order Amounts

Initial Order: Single use bags- 3 packs of 10 (30)

Reusable bags 1 pack of 10

Subsequent Orders: Single use bags – when down to 10 bags reorder 3 packs (30)

Reusable bags – when down to 2 bags reorder 1 pack of 10

Catheter Fixation

Inadequate securement of the catheter and bag can cause pain, bleeding, urethral damage, catheter dislodgement and increase infection risks.

Straps and sleeves are all washable.

Leg straps should be applied tightly enough to prevent the bag slipping down and causing traction on the catheter, but not so tightly that skin damage and circulatory impedance may occur.

Where there is a risk of leg straps being applied incorrectly, the patient should be offered a sleeve-style bag support for which the leg must be measured.

‘All in one’-style sleeves are easier to apply than the alternative style.

The Simpla G-strap is useful for anchoring the inlet tube of long tube leg bags to the thigh. Where peripheral circulation is impaired, an adhesive product such as Clinifix may be more suitable.

Catheter Fixation

Manufacturer	Product	Pack Size	Price (£)	Order Code
Coloplast Ltd	Simpla G-Strap	5	14.68	383001
Bard Ltd	StatLock Catheter Stabilising Device		13.02	FOL01102DT
Wellform Medical	Clinifix Tube Holder: 110mm x 25mm 140mm x 40mm	10	17.69	40-310
		10	19.34	40-410
Optimum Medical Solutions	Ugo fix gentle	5	12.53	3004
Flexicare Medical Ltd	Leg Bag Straps (pairs), washable, anti-slip with velcro fastening	10	11.70	00-0032C

Recommended Order Amounts

Initial Order: 1 box of 10 or 2 boxes of 5

Subsequent Orders: reorder when down to 2

Catheter Bag Sleeve

Manufacturer	Product	Pack Size	Price (£)	Order Code
Coloplast Ltd	Aquasleeve: Small (leg circum. 24-33cm)	4	8.64	783678
	Standard (leg circum. 34-39cm)			
	Medium (leg circum. 40-46cm)			
	Large (leg circum. 47-64cm)			
	Extra Large (leg circum. 65cm+)			

Recommended Order Amounts

Initial Order: 2 packs of 4 (8)

Subsequent Orders: as needed

Catheter Maintenance Solution

The maintenance solution should flow into the bladder by gravity i.e. by elevating the bag above bladder height; it should NOT be squeezed in.

Sodium chloride (NaCl) is used for gentle mechanical flushing of debris or clots.

Suby G solution should only be used for established encrustation and be evaluated for effectiveness and adverse effects on an individual basis.

'Uro-tainer Twin' provides two sequential administrations of a small volume of solution which is more effective than a single larger volume administration.

Uro-Tainer® Polihexanide 0.02% is used for routine decolonization (removal of bacteria) of the catheter by mechanical rinsing.

Catheter Maintenance Solution

Manufacturer	Product	Pack Size	Price (£)	Order Code
Braun Medical Ltd	Uro-tainer NaCl (0.9%)	50ml	3.57	FB99849
		100ml	3.57	FB99833
	Uro-tainer Twin Suby-G (3.23% citric acid)	2x30ml	4.89	9746609
	Uro-tainer Suby-R (6% citric acid)	2x30ml	4.89	9746625
	Uro-tainer PHMB (0.02% Polihexanide)	100ml	3.46	FB99965

Recommended Order Amounts

Initial Order: 7 (but order depending on frequency of use)

Subsequent Orders: dependant on patient need

Catheterisation Trays

These are for Nursing Homes only & should not be ordered for use in the community.

Catheterisation Trays

Manufacturer	Product	Size	Price (£)	Box size	Order Code
Bard Ltd	Bard Bardex I.C. Bard Tray	Male			
		Small 12-16ch	18.32	10	D23655S12-16
		Medium 12-16ch	18.32	10	D23655M12-16
		Large 12-16ch	18.32	10	D23655L12-16
		Female			
		Small 12-16ch	18.32	10	D23695S12-16
		Medium 12-16ch	18.32	10	D23695M12-16
		Large 12-16ch	18.32	10	D23695L12-16
		Bard Biocath Comprehensive Care Foley Tray	Male		
	Small 12-16ch		17.40	10	D22655S12-16
	Medium 12-16ch		17.40	10	D22655M12-16
	Large 12-16ch		17.40	10	D22655L12-16
	Female				
	Small 12-16ch		17.40	10	D22695S12-16
	Bard Lubri-Sil Comprehensive Care Foley Tray	Male			
		Small 12-16ch	18.37	10	D17585S12-16
		Medium 12-16ch	18.37	10	D17585M12-16
		Large 12-16ch	18.37	10	D17585L12-16
Female					
Small 12-16ch		18.37	10	D17615S12-16	
	Medium 12-16ch	18.37	10	D17615M12-16	
	Large 12-16ch	18.37	10	D17615L12-16	

Recommended Amounts

Initial Order: 1 box of 10 per patient

Subsequent Orders: when patient is down to 2 trays

Sheaths

Urinary sheaths are used for males with moderate to severe urinary incontinence or frequency.

Each patient must be assessed for penile girth & length.

Sheath fitting requires dexterity and motivation; success is more likely if the fitting technique is demonstrated.

Sheaths are changed after 24 hours for hygiene purposes and leg/night bags changed every 5 to 7 days.

Correct fitting of the leg bag is important to the success of the sheath and ideally the leg bag should be fitted to the inside of the calf to prevent pooling of urine in the sheath.

Night bags can be connected directly to the sheath; they do not have to be connected to the leg bag, as in the case of catheters.

Urinary Sheaths

Manufacturer	Product	Pack size	Price (£)	Order Code
Coloplast Ltd	Conveen Optima Shorter length:			
	21mm	30	53.07	22121
	25mm	30	53.07	22125
	30mm	30	53.07	22130
Coloplast Ltd	Conveen Optima Standard length:			
	25mm	30	53.07	22025
	28mm	30	52.57	22028
	30mm	30	53.07	22030
	35mm	30	53.07	22035
Bard Ltd	Clear Advantage Spirit:			
	Style 1: Standard			
	24mm	30	51.89	5243
	28mm	30	51.89	5283
	32mm	30	51.89	5323
	36mm	30	51.89	5363
	40mm	30	51.89	5403
	Style 2: Pop-on			
	24mm	30	51.89	7243
	28mm	30	51.89	7283
	32mm	30	51.89	7323
	36mm	30	51.89	7363
	40mm	30	51.89	7403
	Style 3: Wideband			
	24mm	30	51.89	9243
	28mm	30	51.89	9283
32mm	30	51.89	9323	
36mm	30	51.89	9363	
40mm	30	51.89	9403	
Bard Ltd	Clear Advantage with Aloe:			
	Style 1: Standard			
	24mm	30	52.30	1243
	28mm	30	52.30	1283
	32mm	30	52.30	1323
	36mm	30	52.30	1363
	40mm	30	52.30	1403
	Style 2: Pop-on			
	24mm	30	52.30	2243
	28mm	30	52.30	2283
	32mm	30	52.30	2323
	36mm	30	52.30	2363
	40mm	30	52.30	2403
	Style 3: Wideband			
	24mm	30	52.30	3243
	28mm	30	52.30	3283
32mm	30	52.30	3323	
36mm	30	52.30	3363	
40mm	30	52.30	3403	

Recommended Order Amounts**Initial Order:** 1 box of 30**Subsequent Orders:** when down to 10 sheaths reorder 1 box of 30Leg Bags
for Sheaths

Manufacturer	Product	Capacity / Tap	Pack Size	Price (£)	Order Code
Coloplast Ltd	Simpla Profile Leg Bag with anti-kink tubing	500ml / lever lock	10	30.86	21584
Coloplast Ltd	Conveen Active adjustable inlet tube	250ml	10	28.47	25501

Recommended Order Amounts**Initial Order:** 1 box of 10**Subsequent Orders:** reorder 1 box when 2 bags are left

Night Bags for Sheaths

Manufacturer	Product	Capacity / Tap	Pack Size	Price (£)	Order Code
Unomedical Ltd	Careline Easi MT drainage bag with 90cm inlet tube	2L / twist off tube outlet	10	3.40	47-60-LBH
Teleflex	Reusable night bag with integral hanger, standard non-return valve and vent	2L / Slide 2L / Lever	10 10	12.15 12.15	850420-000000 850421-000000

Recommended Order Amounts

Initial Order: Single use bags- 3 packs of 10 (30)
Reusable bags 1 pack of 10

Subsequent Orders: Single use bags – when down to 10 bags reorder 3 packs (30)
Reusable bags – when down to 2 bags reorder 1 pack of 10

Skin Preparation for Sheath Use

Manufacturer	Product	Pack Size	Price (£)	Order Code
Coloplast Ltd	Conveen Prep	54 wipes	10.90	62042

**For wipes & spray adhesive remover request from the company when placing an order*

Recommended Order Amounts

Initial Order: 1 box of 54

Subsequent Orders: reorder when 7 wipes left

Sheath Fixation Devices

Manufacturer	Product	Pack Size	Price (£)	Order Code
Coloplast Ltd	Simpla G-Strap	5	14.68	383001
Bard Ltd	StatLock Catheter Stabilising Device	5	13.02	FOL0102DT
Wellform Medical	Clinifix Tube Holder: 110mm x 25mm 140mm x 40mm	10 10	17.69 19.34	40-310 40-410
Optimum Medical Solutions	Ugo fix gentle	5	12.53	3004
Flexicare Medical Ltd	Leg Bag Straps (pairs), washable, anti-slip with velcro fastening	10	11.70	00-0032C

Recommended Order Amounts

Initial Order: 2 boxes of 5 (10) or 1 box of 10

Subsequent Orders: dependant on patient need

Alternative to Sheaths

Manufacturer	Product	Pack size	Price (£)	Order Code
CliniMed Ltd	BioDerm xls	10	85.75	20026/10
Hollister Ltd	Male retracted penis pouch	10	29.96	9811

Recommended Order Amounts

Initial Order: under advice from the continence advisors

Subsequent Orders: under advice from the continence advisors

Intermittent Catheters

Manufacturer	Product	Length / Size	Price (£)	Box size	Order Code
Teleflex	Liquick Base	Male 10-16ch	85.24	60	630010-630016
		Female 6-14ch	85.24	60	630106-630114
		Tiemann Tip 12-16ch	85.24	60	631012-631016
	Liquick Pure	Female Nelaton tip 8-14ch	45.36	30	851321-000080-000140
		Female Ergothan Tip 8-14ch	46.27	30	851322-000080-000140
	Liquick X-treme	Male (40cm) 12-16ch	50.40	30	851842-000120-000160
	Female (20cm) 12-14ch	50.40	30	851822-000120-000140	
Hollister	VaPro Hydrophilic catheter with integrated hydrating vapour strip	Male (40cm) 8-16ch	53.18	30	72084-72164
		Female (20cm) 8-14ch	53.18	30	72082-72142
		Tiemann Tip 12-16ch	52.71	30	73124-73164
	VaPro pocket Hydrophilic Intermittent Catheter	Male (40cm) 10-16ch	52.68	30	70104-70164
		Female (20cm) 10-14ch	52.82	30	7102-70142
	Infyna Hydrophilic Intermittent catheter	Female (17cm) 8-14ch	47.16	30	88081-88141
		VaPro Plus	Male (40cm) 12-14ch	92.56	30
	Female (20cm) 12-14ch		92.56	30	4122-74142
	Infyna Chic Hydrophilic	Female 8-14ch	50.40	30	7008-7014
Coloplast	SpeediCath Compact Eve	Female 10-14ch	51.30	30	28110-28114
	SpeediCath Compact Plus	Female 10-14ch	49.93	30	28810-28814
	SpeediCath Flex	Male 10-16ch	55.10	30	2892028924
	SpeediCath Compact Set	Female 10ch	53.54	20	28520
		12ch	53.54	20	28522
14ch		53.54	20	28524	
	Male 12-18ch	53.54	20	28422	
Bard	HydroSil gripper	Male 10-18ch	46.05	30	73610-73618
	HydroSil Go	Male 10-20ch	51.50	30	73810G-73820G
		Female 8-18ch	47.67	30	7188-71818

Intermittent Catheters Continued

Manufacturer	Product	Length / Size	Price (£)	Box size	Order Code
ConvaTec	GentleCath Glide Hydrophilic catheter with water sachet	Male 10-16ch	37.89	30	421565-421568
		Female 10-16ch	37.89	30	421571-421578
	GC Glide Hydrophilic catheter with water sachet	Male 8ch,18ch	37.50	30	421564, 421569
		Female 8ch	37.50	30	421570
		Tiemann Tip 8-18ch	37.50	30	421907-421912

Recommended Order Amounts

Initial Order: dependant on frequency of catheterisation

Subsequent Orders: dependant on frequency of catheterisation

Meatal Dilator

Manufacturer	Product	Size	Price (£)	Box size	Order Code
Hunter Urology Ltd	Aloe Meatal	10-16ch	40.30	30	ALOEM10-ALOEM16
	Aloe Meatal Plus	10-16ch	40.51	30	ALOEMPS10-ALOEMPS16

Recommended Order Amounts

Initial Order: dependant on frequency of catheterisation

Subsequent Orders: dependant on frequency of catheterisation

Anal Irrigation

Manufacturer	Product	Contents	Price (£)	Order Code
Coloplast	Peristeen System	1 control unit, 2 rectal catheters, 1 water bag & 2 straps	77.53	29121
	Peristeen Accessory Unit	15 rectal catheters & 1 water bag	135.12	29122
	Peristeen Rectal Catheters	10 rectal catheters	86.39	29123
	Peristeen Straps	Pair	5.54	29124
	Peristeen Tubes	Pack of 2	7.75	29125

Peristeen is the 'Gold Standard' anal irrigation system to be used; recommended in NICE guidelines for managing bowel dysfunction.

Ordering

Any of these companies can be used for all products, however it is recommended to use the companies own home delivery service as stock is guaranteed priority.

Manufacturer	Home Delivery Service	Telephone Number
Bard	Script-easy	0800 470 1930
Coloplast	Charter Healthcare	0800 374 654
Hollister	Fittleworth	0800 378 413
Teleflex	Teleflex Care at Home	0800 032 3383